

Baxter



VIDIAPORT
CARRYING ARM

VIDIAPORT CARRYING ARM

VidiaPort Carrying Arm can be purchased as a stand-alone ceiling mounted, wall mounted and pendant mounted version, or in combination with a Hillrom surgical light.

SPRING ARMS

	VIDIAPORT SPRING ARM 21	VIDIAPORT SPRING ARM 26	VIDIAPORT SPRING ARM 40
Length	920 mm (36.2")	920 mm (36.2")	920 mm (36.2")
Load Capacity	1.5 - 21 kg (3.3 - 46.3 lb)	13 - 26 kg (28.7 - 57.3 lb)	18 - 40 kg (39.7 - 88.2 lb)
Smallest diameter	253 mm2 (10"2)	510 mm2 (20.1"2)	510 mm2 (20.1"2)

MONITOR HOLDER

	SINGLE MONITOR HOLDER		DOUBLE MONITOR HOLDER	
Horizontal rotation	360°	420°	360°	420°
Handle postion	Below or beside the monitor			
Monitor size	1 x 24" – 32" (610 mm x 813 mm)		2 x 24" – 32" (610 mm x 813 mm)	
Max weight per monitor	Determined by load capacity of the respective spring arm			
Interface	Vesa 100			
Monitor holder weight (bottom mounted handles)	5 kg (11 lb)		8 kg (17.6 lb)	
Monitor holder weight (side mounted handles)	5.5 kg (12.1 lb)		9 kg (19.8 lb)	
Monitor holder storage box weight (optional)	1 kg (2.2 lb) per box			

ADAPTER

VIDIAPORT VESA 100 ADAPTER	
Horizontal Rotation	360°
Load Capacity	10 kg (22 lb)
Adapter weight	3 kg (6.6 lb)
Interface	Vesa 75 or Vesa 100

VIDIAPORT CARRYING ARM


COMBINATION WITH OTHER PRODUCTS

	VIDIAPORT PRE- INSTALLATION SET	ILED 7 LIGHT, TL5000 LIGHT, TL3000 LIGHT	TV HD CAMERA	FCS 300 ELECTRO CUBE	MAVIG X-RAY PROTECTION SHIELD INCL. SPRING ARM
VidiaPort Ceiling Single, Duo, Trio	X		X	X	X
VidiaPort Wall				X	X
VidiaPort Pendant			X	X	X
VidiaPort Springarm		X	X	X	

¹Manufacturer material numbers: E-0T25B05, FA102-TR, E-0T54B01-TR, FA402-TR

²Mavig X-Ray Protection Shield and Spring Arm can be combined with duo, trio and quad ceiling installations of iLED7, TL5000, and TL3000

VIDIAPORT CEILING¹ DUO WITH TWO SINGLE MONITOR HOLDERS




Values are in mm.
Subject to changes.

¹Several ceiling tubes and arm lengths available; select 1-4 arms

²Several arm lengths available


VIDIAPORT CARRYING ARM

VIDIAPORT WALL WITH ONE SINGLE MONITOR HOLDER



VIDIAPORT CARRYING ARM

VIDIAPORT PENDANT¹ WITH ONE SINGLE MONITOR HOLDER



Values are in mm.
Subject to changes.

¹Several arm lengths available

For more information please contact your local Baxter representative.

hillrom.eu

The VidiPort carrying arm and iLED 7 Surgical Light are intended to be used by clinicians and medically qualified personnel. This medical device is a regulated healthcare product which, pursuant to such regulation bears a CE mark. Baxter recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Baxter reserves the right to make changes without notice in design, specifications and models. The only warranty Baxter makes is the express written warranty extended on the sale or rental of its products. Baxter Medical Systems GmbH + Co. KG

Baxter, Hillrom, iLED and VidiPort are trademarks of Baxter International, Inc. or its subsidiaries. EMA-CS71-220002 v2 01/2023