

USEFUL DEVICES FOR DAILY PRACTICE

- STOOLS
- MULTI PURPOSE TABLES
- STEPS
- BOLSTER, CUSHION, CUBE
- UPHOLSTERY PROTECTION AND CARE



STOOLS



MULTI PURPOSE
TABLES



STEPS



BOLSTER,
CUSHION,
CUBE

Useful devices for daily practice

Content

Bolster, cushion, cube2

Upholstery protection and care2

Stools

with fixed height3

with height adjustment4

Multi purpose tables6

Steps7



You will find the whole Dewert product range on our website and in our further catalogues. We look forward to hear from you!



Bolster, cushion, cube

Design:

2 different variants for artificial leather eligible for upholstery covering:

- Skai® Pandoria Plus with velvet-smooth surface
- Skai® Toronto EN with innovative staynu-technology



Art. no. 1113-..*
Half roll bolster,
measures (L. x W. x H.):
570 x 200 x 95 mm




Art. no. 1111-..*
Neck roll,
measures: L. 570 x Ø 150 mm



Art. no. 1112-..*
Knee roll,
measures: L. 570 x Ø 200 mm



Art. no. 1211-..*
Wedge shaped bolster,
measures (L. x W. x H.):
400 x 400 x 100 - 10 mm



Art. no. 1015*
Cube, measures (L x W x H):
400 x 300 x 600 mm

Art. no. 1016*
Cube, measures (L x W x H):
500 x 450 x 400 mm



Art. no. 1314-..*
Cushion, velcro fixation,
ticking with polyester filling,
measures (L. x W.): 400 x 400 mm

Upholstery protection and care



Protection mat

Art. no. 1890

made of PVC,
thickness: 0.8 mm,
measures: 700 x 600 mm



Cleaning kit for Skai products

(Art. no. 050-1000)

for careful cleaning of upholstery

Consists of:

- Spray bottle with 200 ml special Skai cleansing agent
- 1 microfibre care cloth 40 x 40 cm

Stool with fixed height

Model 400-01 chromium-plated

Model 400-00 white powder coated (RAL 9010)

Model 400-03 white aluminium powder coated (RAL 9006)

Model 400-04 grey aluminium powder coated (RAL 9007)

- Not movable
- Frame made of square pipe 30 x 30 x 2 mm
- Seat 420 x 420 mm, moulded upholstery, thickness: 65 mm
- Seat height: 480 mm
- Variants for artificial leather eligible for upholstery covering (see below)
- Load capacity: 200 kg
- Weight: 6 kg



Model 404-01 chromium-plated

- Not movable
- 5-leg-stool, made of round pipe, Ø 520 mm, chromium-plated
- Seat plate made of laminated wood, Ø 350 mm
- Seat flat upholstered, thickness: 35 mm
- Seat height: 490 mm
- Variants for artificial leather eligible for upholstery covering (see below)
- Load capacity: 200 kg
- Weight: 4 kg

Design:

- 2 different variants for artificial leather eligible for upholstery covering:
 - Skai® Pandoria Plus with velvet-smooth surface
 - Skai® Toronto EN with innovative staynu-technology



Useful devices for daily practice

Swivel stools with height adjustment

Swivel stools - Models 410-02, 411-02, 412-02, 414-02

- 5-arm polished aluminium base, Ø 560 mm
- **Seat:**
 - Seat height adjustable by gas spring; around going hand ring release under seat for height positioning
 - Seat base made of laminated wood
- **Load capacity:** 130 kg






411-02K, 411-02, 411-02L; in medium height

Design:

- 2 different variants for artificial leather eligible for upholstery covering:
 - Skai® Pandoria Plus with velvet-smooth surface
 - Skai® Toronto EN with innovative staynu-technology



Mobility	Seat upholstery	Flat Ø 360 mm Thickness: 35 mm	Moulded upholstery, Ø 370 mm Thickness: 70 mm
		Model 410-02	Model 411-02
Movable on 5 high-quality universal rollers, Ø 50 mm; for hard and soft floors <ul style="list-style-type: none"> • Cover made of zinc-die-casting, chromium-plated • Rollers according to legal regulations with brackets, only movable under load 		<ul style="list-style-type: none"> • Height adjustment: 490 mm - 680 mm • Weight: 6,5 kg  <p>Variants:</p> <ul style="list-style-type: none"> • 410-02K Height adjustment: 390 mm - 520 mm K • 410-02L Height adjustment: 540 mm - 780 mm L 	<ul style="list-style-type: none"> • Height adjustment: 530 mm - 720 mm • Weight: 8 kg  <p>Variants:</p> <ul style="list-style-type: none"> • 411-02K Height adjustment: 420 mm - 550 mm K • 411-02L Height adjustment: 580 mm - 820 mm L
		Model 412-02 <ul style="list-style-type: none"> • Height adjustment: 460 mm - 650 mm • Weight: 6,5 kg <div style="border: 1px solid black; border-radius: 50%; width: 60px; height: 60px; display: flex; align-items: center; justify-content: center; margin: 20px auto;"> without picture </div>	Model 414-02 <ul style="list-style-type: none"> • Height adjustment: 500 mm - 690 mm • Weight: 8 kg 

Swivel stools with height adjustment

Swivel stool with foot operated height adjustment - Model 417-02

- 5-arm polished aluminium base
- **Mobility:**
 - Movable on 5 high-quality universal rollers, Ø 65 mm; for hard and soft floors
 - Cover made of zinc-die-casting, chromium-plated
 - Rollers according to legal regulations with brackets, only movable under load
- **Seat:**
 - Seat height adjustable by gas spring: 540 mm - 720 mm; Foot operated height adjustment by push plate at the lower frame
 - Seat base made of laminated wood
 - Moulded upholstery, Ø 370 mm, thickness 70 mm
- **Load capacity:** 130 kg
- **Weight:** 9 kg



Swivel Stool - Model 420-02

- 5-arm polished aluminium base
- **Mobility:**
 - Movable on 5 high-quality universal rollers, Ø 65 mm; for hard and soft floors
 - Cover made of zinc-die-casting, chromium-plated
 - Rollers according to legal regulations with brackets, only movable under load
- With going around foot rest
- **Seat:**
 - Seat height adjustable by gas spring: 550 mm - 740 mm; with release hand-lever under seat
 - Seat base made of laminated wood
 - Moulded upholstery, Ø 420 mm, thickness 70 mm
 - Upholstered backrest, manual height adjustable
- **Load capacity:** 130 kg
- **Weight:** 12 kg



Variants:



420-02K
Height adjustment:
440 mm - 570 mm

K



420-02L
Height adjustment:
600 mm - 840 mm

L

Design:

- 2 different variants for artificial leather eligible for upholstery covering:
 - Skai® Pandora Plus with velvet-smooth surface
 - Skai® Toronto EN with innovative staynu-technology



Useful devices for daily practice

Multi purpose tables

Model 325N-00 white powder coated (RAL 9010)

Model 325N-03 white aluminium powder coated (RAL 9006)

Model 325N-04 grey aluminium powder coated (RAL 9007)



Equipment:

- Movable on 4 castors Ø 65 mm
- 2 table-tops made of HPL, thickness 8 mm, white
- Big table top (580 x 400 mm) manual tiltable
- Small table top (150 x 200 mm) fixed
- 1 bag-hook under big table top

Measures:

- Manual height adjustment: 750 - 1.100 mm
Option: Table lifting unforced by gas spring (art. no. 325-55)
- Measures overall: L: 740 mm x W: 400 mm
- Weight: 11 kg



Model 327-00 white powder coated (RAL 9010)

Model 327-03 white aluminium powder coated (RAL 9006)

Model 327-04 grey aluminium powder coated (RAL 9007)

Equipment:

- Movable on 4 castors Ø 65 mm
- 1 table-top, 760 x 420 mm, made of HPL, thickness 8 mm, white
- 1 dustbin, white plastics, removable

Measures:

- Manual height adjustment: 720 - 1.070 mm
Option: Table lifting unforced by gas spring (art. no. 327-55)
- Measures overall: L: 760 mm x W: 420 mm
- Weight: 13 kg

Model 329-00 white powder coated (RAL 9010)

Model 329-03 white aluminium powder coated (RAL 9006)

Model 329-04 grey aluminium powder coated (RAL 9007)

Equipment:

- Movable on 4 castors Ø 65 mm, 2 castors with lock
- 2 table-tops, made of HPL, thickness 8 mm, white
- Upper table-top (measures: 650 x 400 mm) height adjustable: 720 mm - 1.020 mm
- Lower table-top (measures: 580 x 400 mm) fixed

Measures:

- Manual height adjustment of upper table top: 720 mm bis 1.070 mm
Option: Table lifting unforced by gas spring (art. no. 329-55)
- Measures overall: L: 650 mm x W: 400 mm
- Weight: 15 kg



Useful devices for daily practice

Steps with handrail

Model 323-00 white powder coated (RAL 9010)

Model 323-03 white aluminium powder coated (RAL 9006)

Model 323-04 grey aluminium powder coated (RAL 9007)

Equipment:

- With **one-sided handrail** for help and safety of your patients
- Without extra charge: Handrail right or left

Technical data:

- Frame made of square tube 35 x 35 x 3 mm
- Platforms: PUR-Foam with wooden basis; colour: black
- Frame additionally with lateral support to prevent the risk of tipping over



Measures:

- Height step 1: 230 mm | Height step 2: 460 mm | Platforms: 550 x 300 mm
- Height upper edge handrail: 1.350 mm | Width incl. lateral support: 800 mm | Depth: 600 mm

Model 324-00 white powder coated (RAL 9010)

Model 324-03 white aluminium powder coated (RAL 9006)

Model 324-04 grey aluminium powder coated (RAL 9007)

Equipment:

- With **both-sided handrail** for help and safety of your patients
- Without extra charge: Handrail right or left

Technical data:

- Frame made of square tube 35 x 35 x 3 mm
- Platforms: PUR-Foam with wooden basis; colour: black
- Frame additionally with lateral support to prevent the risk of tipping over

without
picture

Measures:

- Height step 1: 230 mm | Height step 2: 460 mm
- Platforms: 550 x 300 mm
- Height upper edge handrail: 1.350 mm
- Width incl. lateral support: 1.050 mm | Depth: 600 mm

Steps made of square tube

Technical data:

- Frame made of square tube, 35 x 35 x 3 mm
- Platforms: PUR-Foam with wooden basis and smooth surface; colour: black
- Platform(s): 550 x 300 mm

*Add colour number for frame:

-01 chromium-plated

-00 white powder coated (RAL 9010)

-03 white aluminium powder coated (RAL 9006)

-04 grey aluminium powder coated (RAL 9007)

Measures **319-..***

Height: 150 mm

Weight: 4,5 kg

Measures **320-..***

Height: 250 mm

Weight: 5 kg

Measures **321-01** (chromium-plated)

Height Step 1: 230 mm | Height Step 2: 460 mm

Length overall: 560 mm | Width overall: 300 mm

Weight: 11 kg



K.H. Dewert GmbH
Vollmestraße 7
33649 Bielefeld
Germany

Post office 14 40 28
33634 Bielefeld
Germany

Phone + 49 (0) 5 21 / 4 00 27-0
Fax + 49 (0) 5 21 / 4 00 27-27

Internet: www.khdewert.de
E-mail: info@khdewert.de
YouTube: K.H. Dewert

05.2021_03

USEFUL DEVICES FOR DAILY PRACTICE

