

CARRYLINE MOBIEL

[OVERZICHT](#)[TECHNISCHE DATA](#)[ACCESSOIRES](#)[DOWNLOADS](#)

APPARATUUR

✓ Verkrijgbaar met vast of in hoogte verstelbaar chassis

- ✓ Synchroonverstelling van het ligvlak (handmatig of elektrisch)
- ✓ De afzonderlijke en handmatige aanpassing van het beenkussen is mogelijk
- ✓ Grote wielen vooraan en een wendbaar chassis voor comfortabel transport
- ✓ 300 mm geleidewielen voor en 125 mm zwenkwielen met vergrendeling achter
- ✓ Comfortabel in- en uitstappen door intelligente functies
- ✓ Automatische verstelling van de rugleuning en beensteun met 10° biedt een ontspannen ligpositie en voorkomt dat de patiënt uit de bank glijdt

POSITIES

The transport chair Carryline Mobil is Available with fixed or height adjustable chassis. The big wheels in the front and an agile chassis ensure a pleasant transport.







MOBILE CHAIR FOR TRANSFER

carryLine mobil in a version with 300mm guide wheels at front and 125mm castors with locking device at rear. The seat height is 51cm.

TRANSPORT CHAIR WITH INDIVIDUAL BENEFITS

The configuration of this transport chair can be perfectly adapted to your specific field of application. The backrest with an anatomical shape and lateral support can be adjusted synchronously with the leg pad - manually or with electromotive assistance. Different backrests can be combined to optimize the transport chair for every individual use.



STABLE STEEL FRAME MADE IN GERMANY

Basis of all carryLine models is a steel frame. CarryLine loungers have a low seat level and a calfrest with negative inclination to allow patients to get in and out comfortably. The automatic adjustment of the seat by 10° offers a relaxed horizontal position, prevents the patient from slipping out of position and relieves the strain of the knees. As the armrests may fully be fold in, free access is ensured.

COVERAGE, COLOUR AND DESIGN

Numerous colour combinations are available for the transport chair, which you can put together individually and adapt to your practice or hospital furnishings. The patient-friendly sitting and lying comfort is always the focus of the upholstery and cover materials. The material of the chair fulfills the requirements for comfort and appearance as well as the daily routine for durability, cleaning and disinfection.



Specifically, GREINER cover material is available in two versions: First, the "compact." Line, which is characterized by patient comfort, longevity and a wealth of variants in 20 colours. And secondly, the innovative "premium." line, which sets new standards with its super-soft and at the same time far more durable surface, and which, in combination with pressure-relieving springs in the seat section, continues to provide unrivaled comfort to this day. Premium cover materials for the transport chair are available in 12 different colours.