

PENTALED 12/28

UNPARALLELED PERFORMANCE TECHNOLOGY,
THE BEST LAMP FOR OUTPATIENT CLINICS AND MINOR SURGERY.

Their compact dimensions and extremely handy structures are permeated by the strong determination of RIMSA's dedication, a synonym of high technology, quality and performance since forever. The thin dome, with its two convenient side grips, ensures easier positioning and adjustment and reduces overall dimensions.

Choose the lamp that best meets your necessities: PentaLED 12 comes with indirect light; PentaLED 28 on the other hand adopts a direct light technology and has the possibility to manually adjust the light field diameter.



PENTALED 12

PENTALED 28



CEILING SINGLE
PENTA12SO
PENTA28SO



MOBILE
PENTA12PI
PENTA28PI



WALL MOUNTED
PENTA12PA
PENTA28PA



CEILING DOUBLE
PENTA12+12
PENTA28+28

Performances


	12	28
Illumination (Ec) at 1 m distance	100.000 lux	120.000 lux
Light head diameter	40 cm	40 cm
Colour temperature	4.500 K	4.500 – 5.000 K
Diameter adjustment	Fixed	Manual
Light field diameter	160 mm	110 – 330 mm
Average LED life	> 60.000 hours	> 60.000 hours
Light intensity control	20 – 100 %	20 – 100 %

Main features

Indirect light – PENTALED 12

Direct light – PENTALED 28

Accessories



Battery group