

# UNICA

THE ONLY LAMP THAT COMPLETELY ELIMINATES GLARE

With the models of the Unica series, RIMSA has obtained the complete elimination of glare. Thanks to its double reflection technology, it is possible to obtain a complete suppression of shadows, the maximization of the light emitting surface and an absolute glare free light.

UNICA 860  
UNICA 520





CEILING SINGLE  
UNICA520SO



MOBILE  
UNICA520PI



CEILING DOUBLE  
UNICA520+520



CEILING SINGLE  
UNICA860SO



CEILING DOUBLE  
UNICA860+520

## Performances

	860	520
Illumination (Ec) at 1 m distance	160.000 lux	160.000 lux
Light head diameter	86 cm	52 cm
Colour temperature (7 selections / continuous)	3.800 to 5.400 K	3.800 to 5.400 K
Light field diameter	170 – 400 mm	170 – 380 mm
Average LED life	> 60.000 hours	> 60.000 hours
Light intensity control	3 – 100 %	3 – 100 %

## Main features

Double reflection: zero glare

Deep 3D light

Compensation area

Light always on focus

## Accessories




Double yoke




Camera:  
FULL-HD cable  
4K cable



Monitor



Battery group



Wall control  
panel



IR remote  
control



Integration  
system